

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human PCID2 protein



## Basic Information

Catalog Number:

Ag24417

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET30a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MRLTDVVQQLVYEAIDSRDGASCAELVSFKHPHVAN  
PRLQMASPEEKCCQVLEPPYDEMFAAHLRCTYAVGD  
HDFIEAYKCQTVIVQSFLRAFQAHKEENWALPVMYA  
VALDLRVFANNADQQLVKKGSKVGDMLEKAAELLM  
SCFRVCASDTRAGIEDSKKWGMLFLVNQLFKIYFKINK  
LHLCKPLIRAISSNLKDDYSTAQRVTYKYVGRKAMF  
DSDFKQAEYLSFAFEHCHRSSQKNRMILYLLPVKM  
LLGHMPTVELKKYHLMQFAEVTRAVSEGNLLLHEAL  
AKHEAFFIRCIFLILEKLIITYRNLFKKVYLLKTHQLS  
LDAFLVALKFMQVEDVDIDEVQCILANLIYMGHVKG  
YISHQHKLKLVVSKQNPFPPLSTVC  
(1-397 aa encoded by BC031246)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : [proteintech@ptglab.com](mailto:proteintech@ptglab.com) W:  
[ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.