

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human PRND protein



Basic Information

Catalog Number:

Ag25523

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

G I K H R I K W N R K A L P S T A Q I T E A Q V A E N R P G A F I K Q G R
K L D I D F G A E G N R Y Y E A N Y W Q F P D G I H Y N G C S E A N V T
K E A F V T G C I N A T Q A A N Q G E F Q K P D N K L H Q Q V L W R L
V Q E L C S L K H C E F W L E R G A G L R V T M H Q P V L L C L L A L I W
L M V K

(28-176 aa encoded by BC043644)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com W:
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.