

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human FLII protein



Basic Information

Catalog Number:

Ag26865

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET30a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

DQAKQVLKGMSDVAQEKNNKQEEADARAPSGKVR
RWDQGLEKPRLDYSEFFTEDVGQLPGLTIWQIENFVP
VLVEEAFHGKFYEADCYIVLKTFLDDSGSLNWEIYYWI
GGEATLDKKACSAIHAVNLRNYLGAECRTVREEMGDE
SEEFQVFDNDISYIEGGTASGFYTVEDTHYVTRMYR
VYGKKNIKLEPVPLKGTSLDPRFVFLDRGLDIYVWRG
AQATLSSTTKARLFAEKINKNERKGKAEITLLVQGQEL
PEFWREALGGEPEIKKHVPEDFWPPQPKLYKVLGLG
YLELPQINYKLSVEHKQRPKVELMPRMRLQSLDTRC
VYILDWSDVFIWLGRIKSPRLVRAAALKLGQELCGML
HRPRHATVSRSLGTEAQVFKAFFKNWDDVLTVDYT
RNAEAVLQSPGLSGKVKRDAEKKDQMKADLTALFLPR
QPPMSLAEEQLMEEW

(440-900 aa encoded by BC025300)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.