

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# HIST3H2BB Fusion Protein



## Basic Information

Catalog Number:

Ag26901

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MPDPSKSA P A P K K G S K K A V T K A Q K K D G K K R K R G R K E S  
Y S I V Y V K V L K Q V H P D T G I S S K A M G I M N S F V N D I F E R I A  
S E A S R L A H Y N K R S T I T S R E V Q T A V R L L L P G E L A K H A V S  
E G T K A V T K Y T S S K  
(1-126 aa encoded by BC100855)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.