

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human CHST1 protein



## Basic Information

Catalog Number:

Ag3183

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

HIILATTRSGSFVQGLFNQHLDFYLFEPYHVQNT  
LI PRFTQ GKSPADRRV MLGASRDLLRSLYDCDLYFLEN  
YIKPPPVNHTTDRI FRRGASRVLC SRPVCDPPGPADLV  
LEEGDCVRKCGLLNLTVAEEACRERSHVAIKTVRVPE  
VNDLRALVEDPRLNLKVIQLVRDPRGILASRSETFRDITY  
RLWRLWYGTGRKPYNLDTVTQTTVCEDFSNSVSTGL  
MRPPWLKGYMLVRYEDLARNPMKKTEEIYGLGIPLD  
SHVARWIQNNTRGDPTLGKHKYGTVRNSAATAEKWR  
FRLSYDIVAFAQNACQVLAQLGYKIAASEEELKNPS  
VSLVEERDFRPFS

(62-411 aa encoded by BC022567)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.