

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

RPL4 Fusion Protein



Basic Information

Catalog Number:

Ag32272

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

HQTSAESWGTGRAVARIPRVGGGTHRSGQGAFGN
MCRGGRMFAPTKTWRRWHRRVNTTQKRYAICSALA
ASALPALVMSKGHRIIEVPELPLVVEDKVEGYKKTKEA
VLLKLLKAWNDIKKVYASQRMRAKGKMRNRRRIQR
RGPCIIYNEDNGIHKAFRNIPGITLLNVSKLNLKLA
PGG HVG R F C I W T E S A F R I K L D E L Y G T W R K A A S L I K S N Y N L P M
H K M I N T D L S R I L K S P E I Q R A L R A P R K K I H R R V L K K N P L K N
L R I M L K L N P Y A K T M R R N T I
(60-340 aa encoded by BC009888)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.