

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

TSP50 Fusion Protein



Basic Information

Catalog Number:

Ag3247

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

LR SAGC W G A G E A P G A L S T A D P A D Q S V Q C V P K A T C P
S S R P R L L W Q T P T T Q T L P S T T M E T Q F P V S E G K V D P Y R S
C G F S Y E Q D P T L R D P E A V A R R W P W M V S V R A N G T H I C
A G T I I A S Q W V L T V A H C L I W R D V I Y S V R V G S P W I D Q M
T Q T A S D V P V L Q V I M H S R Y R A Q R F W S W V G Q A N D I G L
L K L K Q E L K Y S N Y V R P I C L P G T D Y V L K D H S R C T V T G W G L
S K A D G M W P Q F R T I Q E K E V I I L N N K E C D N F Y H N F T K I P
T L V Q I I K S Q M M C A E D T H R E K F C Y E L T G E P L V C S M E G T
W Y L V G L V S W G A G C Q K S E A P P I Y L Q V S S Y Q H W I W D C
L N G Q A L A L P A P S R T L L L A L P L P L S L L A A L
(32-385 aa encoded by BC033016)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.