

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human AD7c-NTP protein



Basic Information

Catalog Number:

Ag32506

Size:

50 µg

Form:

Available lyophilized

Species:

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MEFSLLLRLECNAGISAHRNLRPGSSDPASASPVA
GITGMCTHARLIYFFLVEMEFHVHGQAGLELPTSDDP
SVSASQSARYRTGHHARLCLANFCGRNRVSLMCP
SW SPELKQSTCLSLPKCWDYRRAAVPGLFILFLR
HRCPTL TQDEVQWCDHSSLQPSTPEIKHPPASASQVAGTKD
MHHTWLIFFIFNLRQSLNSVTQAGVQWRNLGSLQ
PLPPGFKLFSCPSLLSSWDYRRPPRLANFFVFLVEMGF
TMFARLILISGPCDLPASASQSAKITGVSHARLIFNF
CLFEMESHVSTQAGVQWPNLGSQPLPPGLKRFSCLS
LPSSWDYGHLPHPANFCIFIRGGVSPYLSGWSQTPD
LR

(1-375 aa encoded by AAC08737.1)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com W:
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.