

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# XPOT Fusion Protein



## Basic Information

Catalog Number:

Ag3256

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MDEQALLGLNPNADSDFRQRALAYFEQLKISPDWQ  
VCAEALAQRITYSDDHVKFFCFQVLEHQVKYKSEL  
TVQQQLIRETLISWLQAQMLNPQPEKTFIRNKAQV  
FALLFVTEYLTKWPKFFDILSVVDLNPRGVDLYLRILM  
AIDSELVDRDVVHTSEEARRNTLIKDTMREQCIPNLVE  
SWYQILQNYQFTNSEVTCQCLEVVGAYVSWIDLSLI  
ANDRFINMLLGHMSIEVLREEACDCLFEVVNKGMDPV  
DKMKLVESLCQVLQSA GFFSIDQEEDVDFLARFSKLV  
NGMGQSLIVSWSKLIKNGDIKNAQEALQAIETKVALM  
LQLLIHEDDDISSNIIGFCYDYLHILKQLTVLSD  
(1-366 aa encoded by BC020569)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.