

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human AKT3 protein



Basic Information

Catalog Number:

Ag32838

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MSDVTIVKEGWVQKRGEYIKNWRPRYFLKTDGSFIG
YKEKPDVLDLPPYPLNFSVAKCQLMKTERPKPNTFIIR
CLQWTTVIERTFHVDTPEREWEAIQAVADRLQR
QEEERMNCSPTSQIDNIGEEEMDASTTHHKRKTMMND
FDYLLKLGKTFGKVLVREKASGKYAMKILKKEVIAK
DEVAHTLTESRVLKNTRHPFLTSKYSFQTKDRLCFVM
EYVNGGELFFHLSRERVFSEDRTRFYGAEIVSALDYH
SGKIVYRDLKLENLMRDKDGHIKITDFGLCKEGITDAAT
MKTFCGTPEYLAPEVLEDNDYGRAVDWWGLGVVVMY
EMMCGRLPFYNQDHEKLFELILMEDIKFPRTLSSDAKS
LLSGLLIKDPNKRLLGGPDDAKEIMRHSFFSGVNWQD
VYDKLVPFPKQVTSSETDRYFDEEFTAQTITITPPE
KCQQSDCGMLGNWKK
(1-465 aa encoded by NM_005465)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.