

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human CXADR protein



## Basic Information

Catalog Number:

Ag34346

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MALLLCFVLLCGVVDFAFARSLITTPPEEMIEKAKGETAYL  
PCKFTLSPEDQGPLDIEWLISPADNQKVDQVILYSGD  
KIYDDYYPDLKGRVHFTSNDLKSGDASINVTNLQLSDI  
GTYQCKVKKAPGVANKKHLVVLVKPSGARYVDGS  
EEIGSDFKIKCEPKESLPLQYEWQKLSDSQKMPTSW  
LAEMTSSVISVKNASSEYSGTYSCTVRNRVGSQDQCLL  
RLNVVPPSNKAGLIAGAIIGTLALALIGLIIFCCRKKRR  
EEKYEKEVHHDIREDVPPPKSRTSTARSYIGSNHSSLG  
SMSPSNMEGYSKTQYNQVPSEDFERTPQSPTLPPAK  
VAAPNLSRMGAIPVMIPAQSKDGSIV\*  
(1-365 aa encoded by BC003684)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : [proteintech@ptglab.com](mailto:proteintech@ptglab.com) W:  
[ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.