

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human IL13RA1 protein



## Basic Information

Catalog Number:

Ag34616

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MEWPARLCGLWALLCAGGGGGGGAAPTETQPPV  
TNLSVSVENLCTVIWTWNPPEGASSNCSLWYFSHFG  
DKQDKKIAPETRRSIEVPLNERICLVGSQCSTNESEK  
PSILVEKCI SPPEGDPE SAVTELQCIWHNLSYMKCSWL  
PGRNTSPDTNYTLYYWHRSLEKIHQCENIFREGQYFG  
CSFDLTKVKDSSEFEQHSVQIMVKDNAGKIKPSFNIVP  
LTSRVKPDPPHIKNSLFHNDLQVWENPQNFISRCL  
FYEVNNSQTEHNVFYVQEA KCENPEFERNVENT  
SCFMVPGVLPD LNTVRIRVKTNKLCYEDDKLWSNWS  
QEMSI GKKRNSTLYITMLLIVPVI VAGAIIVLLYLKRLKI  
IIFPPI PDPGKIFKEMFGDQND DTLHWKKYDIYEKQTK  
EETDSVLLIENLKKASQ  
(1-427 aa encoded by BC009960)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.