

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human ENOSF1 protein



## Basic Information

### Catalog Number:

Ag34665

### Size:

50 µg

### Form:

Available lyophilized

### Species:

human

### Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

### Biological Activity:

Not tested

### Endotoxin Level:

Please contact the lab for more information

### Peptide Sequence:

MVRGRISRLSVRDVRFPTSLGGHGADAMHTDPDYSA  
AYVVIETDAEDGIKGCGITFTLGKGTEVVCAVNALA  
HHVLNKDLKDIVGDFRGFYRQLTSDGQLRWIGPEKGV  
VHLATAAVLNAVWDLWAKQEGKPVWKLVDMDPR  
MLVSCIDFRYITDVLTEEDALEILQKGQIGKKEREKQM  
LAQGYPAYTTSCAWLGYSDDLKQLCAQALKDGWTR  
FKVKVGADLQDDMRRCQIIRDMIGPEKTLMDANQR  
WDVPEAVEWMSKLAFFKPLWIEEPTSPDDILGHATIS  
KALVPLGIGIATGEQCHNRVIFKQLLQAKALQFLQIDS  
CRLGSVNENLSVLLMAKKFEIPVCPHAGGVGLCELVQ  
HLIIFDYISVSASLENRVCEYVDHLHEHFKYPVMIQRA  
SYMPPKDPGYSTEMKEESVKKHQYPDGEVWKKLLPA  
QEN

(1-443 aa encoded by BC001285)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.