

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human PIAS4 protein



Basic Information

Catalog Number:

Ag37383

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

PNIAVKVNHSYCSVPGYYPSNKPGEKRPCRPINLTH
LMYLSSATNRITVTWGNYGKSYSVALYLVRLTSSELL
QRLKTIGVKHPELCKALVKEKRLDPDEIATTGVRVSL
ICPLVKMRLSVPCRAETCAHLQCFDAVFYLMNEKKP
TWMCPVCDKPAPYDQLIIDGLLSKILSECEDADEIYL
VDGSWCPIRAEKERSCSPOGAILVLGPSDANGLLPAP
SVNGSGALGSTGGGPGVSMENGKPGADVVDLTLDL
SSSSEDEEEEEEEDEEGPRPKRRCPFQKGLVPAC
(210-510 aa encoded by BC029874)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com W:
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.