

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human RAI2 protein



## Basic Information

Catalog Number:

Ag3946

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

LPFEAVLQNLFP SQGT LGPPPC QPPPGYAPVPPQPFS  
SPLSPLVPPATLLV PYPVIVPLVPVPIPIPIV PQSSES  
KFSSSF PKPPSSFGLHPFKGTQTPLEKDELKPF DILQPK  
EYFQLSRHTVIKMGSENEALDLSMKSV PWLKAGEVSP  
PIFQEDAPLDLSVAAHRKSEPPPETLYDSGASVDSSG  
HTVMEKLP SGMEISFAPATSHEAPAMMDSHISSDA  
ATEMLSQPNHPSGEVKAENNIEMVGESQAAKVIVSV  
EDAVPTIFCGKIKGLSGVSTKNFSFKREDSVLQGYDIN  
SQGEESMGNAEPLRKPIKNRSIKLKKVNSQEIHMLPIKK  
QRLATFFPRK

(C-term-349aa encoded by BC027937)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.