

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human TRIM11 protein



Basic Information

Catalog Number:

Ag41944

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a-6*his-self-cleavage, with N-terminal no tag.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MELRTVCRVPGLVETLRRFRGDVTLDPDTANPELI
SEDRRSVQRGDLRQALPDSPERFDPGPCVLGQERFTSGR
HYWEVEVGDRTSWALGVCRENNRKEKELSAGNG
FWILVFLGSYYNSSERALPLRDPPIRVGIFLDYEAGHL
SFYSATDGSLLFIFPEIPFSGTLRPLFSPLSSSPTMTIC
RPKGGSGDTLAPQ

(1-202 aa encoded by BC011629)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com W:
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.