

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human MIP3 β protein



Basic Information

Catalog Number:

Ag4200

Size:

50 μ g

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MA LL L A S L L V L W T S P A P T L S G T N D A E D C C L S V T Q K P I

P G Y I V R N F H Y L I K D G C R V P A V V F T T L R G R Q L C A P P D

Q P W V E R I I Q R L Q R T S A K M K R R S S

(1-98 aa encoded by BC027968)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com W:
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.