

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human SYNGR2 protein



Basic Information

Catalog Number:

Ag4380

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6^oHis.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

ASAFFLVVDAYFPQISNATDRKYLVI GDL LFSALWTF L
W F V G F C F L T N Q W A V T N P K D V L V G A D S V R A A I T F S F F
S I F S W G V L A S L A Y Q R Y K A G V D D F I Q N Y V D P T P D P N T
A Y A S Y P G A S V D N Y Q Q P P F T Q N A E T T E G Y Q P P P V Y
(81-229aa encoded by BC000407)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.