

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.



RECOMBINANT HUMAN MND A

Basic Information

Catalog Number:
Ag4476

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
E. coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Validated Application:
Blocking peptide

Peptide Sequence:
MMNEYKKILLKGFELMDDYHFTSIKSLLAYDLGLTTKMQEEYNR
IKITDLMEKFKQGVACLDKLIELAKDMPSLKNLVNLRKEKSKVA
KKIKTQEKAPVKKINQEEVGLAAPAPTARNKLTSEARGRIPVAQK
RKTPNKEKTEAKRNKVSQEQSKPPGPGSASTSAAVDHPPLPQTS
SSTPSNTSFTPNQETQAGRQVDARRNVPQNDPVTWVWKATAP
FKYESPENGGKSTMFHATVASKTQYFHVKVFDINLKEKFKVRKVMIT
ISDYSECKGMVEIKEASSVDFNQNFVFNRIIEANKTPKISQLY
KQASGTMVYGLFMLQKKSVMHKNTIYEIQDNTGSMVWVSGKW
HNIKCEKGDKLRFLCLQRRTVDRKLRKLVCGSHSFIKVIKAKKNKE
GPMNVN
(1-407 aa encoded by BC032319)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.