

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

ATP2A3 Fusion Protein



Basic Information

Catalog Number:

Ag4527

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

TRRMARKNAIVRSLPSVETLGCTSVICSDKTGTLTTNQ
MSVCRMVVAEADAGSCLLHEFTISGTTYPEGEVR
QGDDQPVRCGQFDGLVELATICALNDSALDYNEAKG
VYEKVGAEATETALTCLVEKMNVDLQALSRVERAG
ACNTVIKQLMRKEFTLEFSRDRKSMSVYCTPTRPHPT
GQGSKMFVKGAPESVIERCSSVRVGSRTAPLTPTSRE
QILAKIRDWGSGSDTLRCLALATRDAPPRKEDMELDD
CSKFVQYETDLTFVGCVGMLDPPRPEVAACITRCYQ
AGIRVVMITGDNKGTAVAICRRLGIFGDTEVAGKAY
TGREFFDLSPEQQRQACRTARCFARVEPAHKSRIVE
NLQSFNEITAMTGDGVNDAPALKKAEIGIAMGSGTA
VAKSAAEMVLSDDNFASIVAAVEEGRAIYSNMKQFIR
YLISSNVGEVVCIFLTAIILGLPEALIPVQLLWVNLVTDG
LPATALGFNPPDLIMEKLPRSPREALISGWLFFRYLAI
GVYVGLATVAAAATWWFVYDAEGPHINFYQLRNFLKC
SEDNPLFAGIDCEVFESRFPTTMALSVLVTIEMCNALN
SVSENSQLLRMPWMPNPLLVAVAMSMALHFLILLV
PPLPLIFQVTPLSGRQWVVV/LQISLPVILLDEALKYLSR
NHMHEEMSQK

(323-999 aa encoded by BC035729)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.