

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human RUFY1 protein



## Basic Information

Catalog Number:

Ag4606

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

IAMKLEKDTHEKQDTLVALRQQLEEVKAINLQMFHK  
AQNAESSLQKNEAITSFEGKTNQVMSSMKQMEERL  
QHSEARQGAEEERSHLQQLGGRIGALQLLSQLH  
EQCSSLEKELKSEKEQRQALQRELQHEKDTSSLLRMEL  
QQVEGLKKEKRELQDEKAEQKICEEQEQALQEMGL  
HLSQSKLKMEDIKEVNQALKGHAWLKDDEATHCRQC  
EKEFSISRKHHCRCNGHIFCNTCSSNELALPSYPKPV  
RVCDSCHTLLLQRCSTAS  
(325-708aa encoded by BC032571)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com) W:  
[ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.