

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human GRM3 protein



## Basic Information

Catalog Number:

Ag4794

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

HNFLRREIKIEGDLVLGGFLPINEKGTGTEECGRINEDR  
GIQRLEAMLFIDEINKDDYLLPGVKLGVHILDTCSRDT  
YALEQSLEFVRASLTQVDEAEYMC PDGSYAIQENIPLL  
IAGVIGGSYSSVSIQVANLLRLFQIPQISYASTSAKLSD  
KSRDYDFARTVPPDFYQAKAMAEILRFFNWTYVSTVA  
SEG DYGETGIEAFEQEARLRNICIATAEKVGRSNIRKS  
YDSVIRELLQKPNARVVVLFMRSDDSRELI AAASRAN  
ASFTWVASD GWGAQESI IKGSEHVAYGAI TLELASQ  
PVRQFD RYFQSLNPYNNHRNPWFRDFWEQKFQCSL  
QNKRNHRRVCDKHLAIDSSNYEQESKIMFVVNA  
(26-396aa encoded by BC022496)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.