

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

VSIG1 Fusion Protein



Basic Information

Catalog Number:

Ag5314

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

FARNKAKAKAKERNKTI A E L E P M T K I N P R G E S E A M P R
E D A T Q L E V T L P S S I H E T G P D T I Q E P D Y E P K P T Q E P A P E
P A P G S E P M A V P D L D I E L E L E P E T Q S E L E P E P E P E S E P
G V V V E P L S E D E K G V V K A
(256-387 aa encoded by BC043216)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com W:
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.