

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human NINJ2 protein



## Basic Information

Catalog Number:

Ag5343

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6<sup>+</sup>His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MESARENIDLPQSSDPRSQPINLNHYATKKSVAESM  
LDVALFMSNAMRLKAVLEQGPSSHYTTLVTLISLSLL  
QVVIGVLLVVIARLNLNEVEKQWRLNQLNNAATILVF  
FTVVINVFITAFGAHKTGFLAARASRNPL  
(17-158 aa encoded by BC057766)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : [proteintech@ptglab.com](mailto:proteintech@ptglab.com) W:  
[ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.