

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human β -ENaC protein



Basic Information

Catalog Number:

Ag5348

Size:

50 μ g

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

GI F I R T Y L S W E V S V S L S V G F K T M D F P A V T I C N A S P F K Y
S K I K H L L K D L D E L M E A V L E R I L A P E L S H A N A T R N L N F S I
W N H T P L V L I D E R N P H H P M V L D L F G D N H N G L T S S S A S E
K I C N A H G C K M A M R L C S L N R T Q C T F R N F T S A T Q A L T E
W Y I L Q A T N I F A Q V P Q Q E L V E M S Y P G E Q M I L A C L F G A
E P C N Y R N F T S I F Y P H Y G N C Y I F N W G M T E K A L P S A N P G
T E F G L K L I D I G Q E D Y V P F L A S T A G V R L M L H E Q R S Y P F I
R D E G I Y A M S G T E T S I G V L V D K L Q R M G E P Y S P C T V N G S
E V P V Q N F Y S D Y N T T Y S I Q A C L R S C F Q D H M I R N C N C
(68-401 aa encoded by BC036352)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com W:
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.