

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human $\beta$ -ENaC protein



## Basic Information

Catalog Number:

Ag5515

Size:

50  $\mu$ g

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

GI F I R T Y L S W E V S V S L S V G F K T M D F P A V T I C N A S P F K Y  
S K I K H L L K D L D E L M E A V L E R I L A P E L S H A N A T R N L N F S I  
W N H T P L V L I D E R N P H P M V L D L F G D N H N G L T S S S A S E  
K I C N A H G C K M A M R L C S L N R T Q C T F R N F T S A T Q A L T E  
W Y I L Q A T N I F A Q V P Q Q E L V E M S Y P G E Q M I L A C L F G A  
E P C N Y R N F T S I F Y P H Y G N C Y I F N W G M T E K A L P S A N P G  
T E F G L K L I D I G Q E D Y V P F L A S T A G V R L M L H E Q R S Y P F I  
R D E G I Y A M S G T E T S I G V L V D K L Q R M G E P Y S P C T V N G S  
E V P V Q N F Y S D Y N T T Y S I Q A C L R S C F Q D H M I R N C N C  
(68-401 aa encoded by BC036352)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.