

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human PILRB protein



Basic Information

Catalog Number:

Ag5909

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

QPGGSTGSGPSYLYGVTQPKHLSASMGGSV EIPFSFY
YPWELAI V PNV R I S W R R G H F H G Q S F Y S T R P P S I H K D Y
V N R L F L N W T E G Q E S G F L R I S N L R K E D Q S V Y F C R V E L D
T R R S G R Q Q L Q S I K G T K L T I T Q A V T T T T T W R P S T T T I A
G L R V T E S K G H S E S W H L S L D T

(20-188aa encoded by BC050547)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.