

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human ASXL2 protein



Basic Information

Catalog Number:

Ag6272

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

APSPTGFQLEDISTSQRFMLGFAGRRTSKPAMAGHYL
LNISTYGRGSESFRRTHSVNPEDRFCLSSPTEALKMGY
TDCKNATGESSSSKEDDTDEESTGDEQESVTVKEEP
QVSQSAGKGDTS SSGPHSRETLSTSDCLASKNVKAEIP
LNEQTTLSKENYLFTRGQTFDEKTLARDLIQAAQKQM
AHAVRGKAI R S S P E L F S S T V L P L P A D S P T H Q P L L L P P L Q
TPKLYGSPTQIGPSYRGMINVSTSSDMDHNSAVPGSQ
VSSNVGDVMSFSVTVTIPASQAMNPSSHGQTIPVQ
AFSEENSI EGT P S K C Y C R L K A M I M C K G C G A F C H D D C I
GPSKLCV S C L V V R

(572-918aa encoded by BC042999)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.