

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human NAP1L1 protein



Basic Information

Catalog Number:

Ag6690

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MADIDNKEQSELDQDLDDVEVEVEEETGEETKLKARQ
LTVQMMQNPQILAAQERLDGLVETPTGYIESLPRVV
KRRVNALIKNLQVKCAQIEAKFYEEVHDLERKYAVLYQ
PLFDKRFEIINAIYEPTEECEWKPDDEEISEELKEKAK
IEDEKKDEEKEDPKGIPEFWLTVFKNVDLLSDMVQEH
DEPILKHLKDIKVKFSDAGQPMSEVLEFHFEPNEYFTN
EVLTKTYMRSEPDSDPF5FDGPEIMGCTGCQIDWK
KGKNVTLTKIKKKQKHKGRTVVRTVTKTVSNDSEFFNF
FAPPEVPESGDLDDAEAILAADFEIGHFLRERIIPRSV
LYFTGEAIEDDDDDYDEEGEEADEEGEEEGDEENDPD
YDPKKDQNPAECKQQ

(1-391 aa encoded by BC002387)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.