

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human NARS protein



## Basic Information

Catalog Number:

Ag6777

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6<sup>o</sup>His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

ELSCDFWELIGLAPAGGADNLNIEESDQVQLNRRH  
MMIRGENMSKILKARSMVTRCFRDHFFDRGYEVTTP  
TLVQTQVEGGATLFKLDYFGEEAFLTQSSQLYLETCLP  
ALGDVFCIAQSYRAEQSRTRRHLAEYTHVEAECPLT  
FDDLNRLEDLVCDDVDRILKSPAGSIVHELNPNFQPP  
KRPFKRMNYSDAIVWLKEHDVKKEDGTFYEFGEDIP  
APERLMTDTINEPILLRFPVEIKSFYMQRCPEDSRLTE  
SVDVLMPNVGEIVGGSMRIFDSEELAGYKREGIDPTP  
YYWYTDQRKYGTCPHGGYGLGLERFLTWILNRYHIRD  
VCLYPRFVQRCTP

(199-548 aa encoded by BC001687)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.