

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human CPT2 protein



Basic Information

Catalog Number:

Ag7071

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

SGNEESLRKVD S A V F C L C L D D F P I K D L V H L S H N M L H G
D G T N R W F D K S F N L I I A K D G S T A I H F E H S W G D G V A V L R
F F N E V F K D S T Q T P A V T P Q S Q P A T T D S T V T V Q K L N F E
L T D A L K T G I T A A K E K F D A T M K T L T I D C V Q F Q R G G K E F L
K K Q K L S P D A V A Q L A F Q M A F L R Q Y G Q T V A T Y E S C S T A
A F I K H G R T E T I R P A S V Y T K R C S E A F V R E P S R H S A G E L Q
Q M M V E C S K Y H G Q L T K E A A M G Q G F D R H L F A L R H L A A
A K G I I L P E L Y L D P A Y G Q I N H N V L S T S L S P A V N L G G F A
P V V S D G F G V G Y A V H D N W I G N C N V S S Y P G R N A R E F L Q
C V E K A L E D V F D A L E G K S I K S

(309-658aa encoded by BC002445)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.