

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human USP10 protein



## Basic Information

**Catalog Number:**  
Ag7106

**Size:**  
50 µg

**Form:**  
Available lyophilized

**Species:**  
human

**Expression Source:**  
*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**  
Not tested

**Endotoxin Level:**  
Please contact the lab for more information

**Peptide Sequence:**  
KFIPLYSKVQRPCSTPMIDSFVRLMNEFTNMPVPPKP  
RQALGDKIVRDIRPGAAFEPTYIYRLLTVNKSSLSEKGR  
QEDAEYLGFI LNGLHEEMLNLKLLSPSNEKLTISNGP  
KNHVSNEEEQEEQEGEGSEDEWEQVGPRNKTSVTRQ  
ADFVQTPITGIFGGHIRSVVYQQSSKESATLQPFRTL  
QLDIQSDKIRTVQDALESLVARESVMQGYTTKTKQVE  
ISRRVTLEKLPVVLVHLKRFVYEKGGCQKLIKNIYYPV  
DLEISKELSPGVKNKFKCHRTRYLFAVVYHHGNSAT  
GGHYTTDFVQIGLNGWLRIDDQTVKVINQYQVVKPT  
AERTAYLLYRRVDLL  
(444-798 aa encoded by BC000263)

## Selected Validation Data

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) E: proteintech@ptglab.com W:  
(toll free in USA), or 1(312) 455-8498 ptglab.com  
(outside USA)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.