

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human TEKT2 protein



## Basic Information

Catalog Number:

Ag7311

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e. coli*-derived, PET28a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MATLSVKPSRRFQLPDWHTNSYLLSTNAQLQRDASH  
QIRQEARVLRNETNNQTIWDEHDNRTRLVERIDTVNR  
WKEMLDKCLDLDAEIDALTMKESAEQNLQAKNLPL  
DVAIECLTLRESRRDIDVVKDPVEDELHKEVEIEATK  
KALQKVSQAFEQLCLLQEVQQQLNSDHRGKMETL  
EIDRGCLSLNLRSPNISLKVDPTRVPDGGSTTLQWDDF  
SRFNKDRAEAEMKAAATELREATAITIAETNNELEAQR  
VATEFAFRKRLREMEKVYSELKWQEKNTLEEIAELQED  
IRHLEEDLRTKLLSLKSHRLEARTYRPNVELCRDQAAQ  
YGLTDEVHQLEATIAALKQKLAQAQDALDALCKHLAR  
LQADIAACKANSMLLDTKCMDTRRRLKTPAERFVPEVD  
TFTRTTNSTLSPLKSCQLELA  
(1-430aa encoded by BC035620)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.