

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human PTPN7 protein



## Basic Information

Catalog Number:

Ag7469

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e. coli*-derived, PET28a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MGASFWPIRQAREQRRALSFRQTSWLSEPLGPAP  
HLSMVQAHGGRSRAQPLTSLGAAMTQPPPEKTPAK  
KHVRLQERRGSNVALMLDVRSLGAVEPICSVNTPREV  
TLHFLRTAGHPLTRWALQRQPPSPKQLEEEFLKIPSNF  
VSPEDLDIPGHASKDRYKTLIPNPQSRVCLGRAQSQE  
DGDYINANYIRGYDGKEKVYIATQGMPMNTVSDFWE  
MVWQEEVSLIVMLTQLREGKEKVHYWPTEEETYPG  
FQIRIQDMKECPEYTVRQLTIQYQEERRSVKHILFSA  
WPDHQPESAGPLLRVAEVEESPETAHGPPIVVHC  
SAGIGRTGCFIATRIGCQQLKARGEVDILGIVCQLRLD  
RGGMIQTAEQYQFLHHTLALYAGQLPEEPS  
(1-399 aa encoded by BC001746)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.