

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human GLA protein



Basic Information

Catalog Number:

Ag7505

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6^{His}.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

WKDAGYEYLCIDDCWMA PQRDSEGR LQADPQRFP H
GIRQLANYVHSKGLKGIYADVGNKTCAGFP GSGFGYY
DIDAQTFADWGV D L L K F D G C Y C D S L E N L A D G Y K H M S
L A L N R T G R S I V Y S C E W P L Y M W P F Q K P N Y T E I R Q Y C N H
W R N F A D I D D S W K S I K S I L D W T S F N Q E R I V D V A G P G G
W N D P D M L V I G N F G L S W N Q Q V T Q M A L W A I M A A P L F
M S N D L R H I S P Q A K A L L Q D K D V I A I N Q D P L G K Q G Y Q L R
Q G D N F E V W E R P L S G L A W A V A M I N R Q E I G G P R S Y T I A
V A S L G K G V A C N P A C F I T Q L L P V K R K L G F Y E W T S R L R S
H I N P T G T V L L Q L E N T M Q M S L
(81-425 aa encoded by BC002689)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com W:
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.