

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human CAB1 protein



## Basic Information

**Catalog Number:**

Ag7509

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6<sup>o</sup>His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Peptide Sequence:**

VRQSADF MPLKQMMKTLNNDLGP NWRDKLEYFEERP  
FAAASIGQVHLARMKGGREVA MKIQYPGVAQ SINS D  
VNNLMAVLNMSNMLPEGLFPEHLIDVLRRELALECDY  
QREAA CARKFRDLLKGHPFFYVPEIVDELCSPHVLTTE  
LVSGFPLDQAEGLSQEIRNEICYNIIVLCLRELF EFHFM  
QTDPNWSNFFYDPOQH KVALLD FGATREYDRSFTDL  
YIQIIRAAADRRET VRAKSIEMKFLTGYEVKVMEDA H  
LDAIILG EAFASDEPFDGTQSTTEKIHNLI PVMLRHR  
LVPPP EETYS LHRKMGG SFLI CSK LKARFPCKAMFEEA  
YSNYCKRQAQQ

(300-647 aa encoded by BC005171)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
[ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.