

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human GPR56 protein



Basic Information

Catalog Number:

Ag7799

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6[°]His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

GHREDFRFC^SQRNQT^HRSSLHYKPT^DLRISIENSEEA
LTVHAPFPAAHPASRSF^PDPGRGLYHFC^LYWNRHAGRL
HLLYGK^RDFLLSDKASSLLCF^QHQEESLAQG^PPLLATS
VTSWWSPQ^NISLPSAASF^TFSFHSPHTAAHNASVD
MCELK^RDLQLLSQ^FLKHPQKASRRPSA^APASQ^LQLQSL
ESK^LTSVRFMGDMVS^FEEDRINATV^WKLQPTAGLQD
LHIHSRQ^EEEQSEIMEYSV^LLPRTLFQ^RTKGRSGEAEK
RLLLVDFSSQ^ALFQDKNSSHV^LGKVLGIVVQ^NTKVA
NLTEPV^LTFQHQLQ^PKNVTLQ^CVFWVEDPTLSSPG
(27-380aa encoded by BC008770)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.