

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human cTnT protein



Basic Information

Catalog Number:
Ag7937

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Peptide Sequence:
MSDIEEVVEEYEEEEEQEEAAVEEQEEAAEEDAEAEA
ETEETRAEEDEEEEAKEAEDGPMEESEKPKPRSFMPNL
VPPKIPDGERVDFDDIHRKRMEKDLNELQALIEAHFEN
RKKEEEELVSLKDRIERRRAERAEQQRIRNEREKERQN
RLAEERARRREEENRRKAEDEARKKKALSMMHFGGY
IQKTERKSGKRQTEREIKKKILAERRKVLAIIDLNLNEDQL
REKAKELWQSIYNLEAEKFDLQEKFKQKYEINVLNRN
INDNQKVSCTRGGKAKVTGRWK
(1-285aa encoded by BC002653)

Selected Validation Data

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.