

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# RECOMBINANT HUMAN GS28

## Basic Information

Catalog Number:

Ag9021

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*E. coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MAAGTSSYWEDLRKQARQLENELDLKLVFSKLCYSYSHSSTR  
DGRFRDRYSDDTTPLLNGSSCDRMFETMAIEEQLLARLTGVNDK  
MAEYTNAGVPSLNAALMHTLQRHRDILQV  
(1-117 aa encoded by BC012620)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA),  
or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.