

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

COPG2 Fusion Protein



Basic Information

Catalog Number:

Ag9058

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6^{His}.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

EVRDRATFYLNVLQQRQMALNATYIFNGLTVSVPGM
EKALHQYTLPESEKPFDMKSIPLAMAPVFEQKAEITLV
ATKPEKLAPSRQDIFQEQLAAIPEFLNIGPLFKSSEPV
QLTEAETEYFVRCIKHMFNHNIVFQFDCTNTLNDQLL
EKVTVQMEPSDSYEVLSVIPAPSLPYNQPGICYTLVRL
PDDDPTAVAGSFSCMKTFRDCDPNTGVPDEEDGY
DDEYVLEDELVTVSDHIQKVLKPNFAAAWEEVGDTF
EKEETFALSSTKTLLEAVNNIITFLGMQPCERSDKVPE
NKNSHSLYLAGIFRGGYDLLVRSRLALADGVMTQVTV
RSKERTPVDVILASVG

(C-term-349aa encoded by BC017443)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.