

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

TCTN3 Fusion Protein



Basic Information

Catalog Number:
Ag9106

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Peptide Sequence:
QTQSPPSFYRAGDPILTYFPKWSVISLLRQPAGVGAG
GLCAESNPAGFLESKSTTCRFFKKNLASSCTLDSALNA
ASYYNFTVLKVPMSMDPQNMEFQVPVILTSQANAP
LLAGNTCQNVVSQVYETNGTFGIQKVSVSLGQT
NLTVEPGASLQQHFILRFRAFQQSTAASLTSRSGNP
GYVIGKPLLALDDISYSMTLLQSQGNGSCSVKRHEV
QFGVNAISGCKLRLLKADC SHLQQEIYQTLHGRPRPE
YVAIFGNADPAQKGGWTRILNRHCSISAINCTSCCLIP
VSLEIQVLWAYVGLLSNPQAHVSGVRFYQCQSIQD
SQQVTEVSLTTLVNFVDITQKPPRGPQPKMDWKW
PFDFFPFKVAFSRGVFSQKCSVSPIL
(1-393 aa encoded by BC009494)

Selected Validation Data

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) E: proteintech@ptglab.com W:
(toll free in USA), or 1(312) 455-8498 ptglab.com
(outside USA)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.