

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

PHKA2 Fusion Protein



Basic Information

Catalog Number:

Ag9495

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

SI AVL TQ EIV VY LAM YV R A Q P S L F V E M L R L R I G L I Q V
M A T E L A R S L N C S G E E A S E S L M N L S P F D M K N L L H H I L S G
K E F G V E R S V R P I H S S T S S P T I S I H E V G H T G V T K T E R S G I
N R L R S E M K Q M T R R F S A D E Q F F S V G Q A A S S S A H S S K S
A R S S T P S S P T G T S S S D S G G H H I G W G E R Q G Q W L R R R R
L D G A I N R V P V G F Y Q R V W K I L Q K C H G L S I D G Y V L P S S T
T R E M T P H E I K F A V H V E S V L N R V P Q P E Y R Q L L V E A I M V
L T L L S D T E M T S I G G I I H V D Q I V Q M A S Q L F L Q D Q V S I G
A M D T L E K D Q A T G I C H F F Y D S A P S G A Y G T M T Y L T R A V
A S Y L Q E L L P N S G C Q M Q

(886-1235aa encoded by BC014036)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com W:
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.