

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

CPM Fusion Protein



Basic Information

Catalog Number:

Ag9761

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

YLVTS D G K D P E I T N L I N S T R I H I M P S M N P D G F E A V K K P
D C Y Y S I G R E N Y N Q Y D L N R N F P D A F E Y N N V S R Q P E T V
A V M K W L K T E T F V L S A N L H G G A L V A S Y P F D N G V Q A T G
A L Y S R S L T P D D D V F Q Y L A H T Y A S R N P N M K K G D E C K N K
M N F P N G V T N G Y S W Y P L Q G G M Q D Y N Y I W A Q C F E I T L
E L S C C K Y P R E E K L P S F W N N N K A S L I E Y I K Q V H L G V K G Q
V F D Q N G N P L P N V I V E V Q D R K H I C P Y R T N K Y G E Y Y L L L L
P G S Y I I N V T V P G H D P H I T K V I I P E K S Q N F S A L K K D I L L P F
Q G Q L D S I P V S N P S C P M I P L Y R N L P D H S A A T
(99-426 aa encoded by BC022276)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com W:
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.