

For Research Use Only

Human IgG Polyclonal antibody

Catalog Number: B900640

1 Publications



Basic Information

Catalog Number:
B900640

Size:
1mg, Concentration: 1 mg/ml by
Bradford method using BSA as the
standard;

Source:
N/A

Isotype:
IgG

GenBank Accession Number:
GeneID (NCBI):

Full Name:
Calculated MW:
155 kDa

Purification Method:
Protein A

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

This product is Normal human IgG which purified with Protein A. The concentration is 1.8ug/ul. An isotype control antibody should have the same immunoglobulin type and be used at the same concentration as the test antibody.

Notable Publications

Author	Pubmed ID	Journal	Application
Junlan Zhong	35200359	Biosensors (Basel)	

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data