

For Research Use Only

Biotin Anti-Human STRO-1 (STRO-1) Mouse IgG2a Recombinant Antibody

Catalog Number: Biotin-65585



Basic Information

Catalog Number:

Biotin-65585

GenBank Accession Number:

GeneID (NCBI):

Purification Method:

Protein A purification

Size:

100ug , 500 ug/ml

Full Name:

CloneNo.:

STRO-1

Source:

Mouse

Isotype:

IgG2a

Applications

Tested Applications:

FC

Species Specificity:

human

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1×10^6 MG-63 cells were surface stained with 1 μ g Biotin Anti-Human STRO-1 (STRO-1) Mouse IgG2a RecAb (Biotin-65585, Clone: STRO-1) and PE-conjugated streptavidin (red), or unstained (blue). Cells were not fixed.