

For Research Use Only

FcZero-rAb™ Biotin Anti-Mouse CD200 Rabbit Recombinant Antibody

Catalog Number: Biotin-FcA98039



Basic Information

Catalog Number:

Biotin-FcA98039

Size:

100ug , 500 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0920

GenBank Accession Number:

/

GeneID (NCBI):

17470

UNIPROT ID:

O54901

Full Name:

CD200 antigen

Calculated MW:

31 kDa

Purification Method:

Protein A purification

CloneNo.:

240391A6

Recommended Dilutions:

FC: 0.25 ug per 10⁶ cells in a 100 µl suspension

Excitation/Emission maxima wavelengths:

-

Applications

Tested Applications:

FC, ELISA

Species Specificity:

mouse

Positive Controls:

FC : mouse splenocytes,

Background Information

CD200, also known as OX2, is a type I transmembrane glycoprotein that belongs to the immunoglobulin superfamily (PMID: 3032785). It contains two extracellular immunoglobulin domains, a transmembrane and a cytoplasmic domain. CD200 is expressed on a variety of cell types, including thymocytes, B lymphocytes, a subset of T lymphocytes, follicular dendritic cells, neurons, and endothelial cells (PMID: 3032785; 10981966). CD200 binds to its receptor, CD200R, which is primarily expressed by myeloid and T cell lineage (PMID: 15187158). The CD200-CD200R interaction plays a role in immunosuppression and suppression of anti-tumor immune responses (PMID: 22020332).

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ mouse splenocytes were surface stained with 0.25 ug Biotin Anti-Mouse CD200 Rabbit RecAb (Biotin-FcA98039, Clone: 240391A6) or 0.25 ug FcZero-rAb™ Biotin Rabbit IgG Isotype Control Recombinant Antibody (Biotin-FcA98136, Clone: 240953C9), and PE Streptavidin (Cat.No. PF00021). Cells were co-stained with APC Anti-Mouse CD45R. Cells were not fixed.



1x10⁶ mouse splenocytes were surface stained with 0.25 ug Biotin Anti-Mouse CD200 Rabbit RecAb (Biotin-FcA98039, Clone: 240391A6) (red) or 0.25 ug FcZero-rAb™ Biotin Rabbit IgG Isotype Control Recombinant Antibody (Biotin-FcA98136, Clone: 240953C9) (blue), and PE Streptavidin (Cat.No. PF00021). Cells were not fixed.

