

For Research Use Only

FcZero-rAb® Biotin Anti-Human CEACAM3/6 (CD66c/d) Rabbit Recombinant Antibody

Catalog Number: Biotin-FcA98169



Basic Information

Catalog Number:

Biotin-FcA98169

Size:

100ug, 500 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0144

GenBank Accession Number:

BC106728

GeneID (NCBI):

1084

UNIPROT ID:

P40198

Full Name:

carcinoembryonic antigen-related
cell adhesion molecule 3

Calculated MW:

252 aa, 27 kDa

Purification Method:

Protein A purification

CloneNo.:

241254D1

Recommended Dilutions:

FC: 0.25 ug per 10⁶ cells in a 100 µl
suspension

Excitation/Emission maxima
wavelengths:

-

Applications

Tested Applications:

FC, ELISA

Species Specificity:

human

Positive Controls:

FC : human peripheral blood leukocytes,

Background Information

Carcinoembryonic antigen-related cell adhesion molecules (CEACAMs) belong to the immunoglobulin superfamily of cell adhesion molecules. CEACAMs play major roles in cell-cell and cell-ECM adhesion and have been implicated in controlling cell proliferation, angiogenesis, and tissue remodeling (PMID: 23021083). CEACAM3, also named as CD66d and CGM1, is exclusively expressed on granulocytes and acts as the major granulocyte receptor mediating recognition and efficient opsonin-independent phagocytosis of CEACAM-binding microorganisms (PMID: 18606569). CEACAM6, also known as CD66c, is normally expressed on the surface of myeloid and epithelial surfaces and upregulated preferentially at the apical/luminal membranes of many tumors (PMID: 27595061; 35082925). This antibody is raised against 35-155 aa of human CEACAM3/CD66d and can cross-react with CEACAM6/CD66c.

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ human peripheral blood leukocytes were surface stained with 0.25 ug Biotin Anti-Human CEACAM3/6 (CD66c/d) Rabbit RecAb (Biotin-FcA98169, Clone: 241254D1) or FcZero-rAb® Biotin Rabbit IgG Isotype Control Recombinant Antibody (Biotin-FcA98136, Clone: 240953C9), and Streptavidin-PE Conjugate (PE-PF00030). Cells were not fixed.



1x10⁶ human peripheral blood leukocytes were surface stained with 0.25 ug Biotin Anti-Human CEACAM3/6 (CD66c/d) Rabbit RecAb (Biotin-FcA98169, Clone: 241254D1) (red) or FcZero-rAb® Biotin Rabbit IgG Isotype Control Recombinant Antibody (Biotin-FcA98136, Clone: 240953C9) (blue), and Streptavidin-PE Conjugate (PE-PF00030). Cells were not fixed.