

For Research Use Only

FcZero-rAb® Biotin Anti-Human DC-SIGNR/CD299 Rabbit Recombinant Antibody

Catalog Number: Biotin-FcA98205



Basic Information

Catalog Number:

Biotin-FcA98205

Size:

100ug, 500 ug/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC038851

GeneID (NCBI):

10332

UNIPROT ID:

Q9H2X3

Full Name:

C-type lectin domain family 4, member M

Calculated MW:

45 kDa

Purification Method:

Protein A purification

CloneNo.:

241683G8

Recommended Dilutions:

FC: 0.25 ug per 10⁶ cells in a 100 µl suspension

Excitation/Emission maxima wavelengths:

-

Applications

Tested Applications:

FC, ELISA

Species Specificity:

human

Positive Controls:

FC : Transfected HEK-293T cells,

Background Information

Dendritic cell-specific intracellular adhesion molecule-3-grabbing non-integrin (DC-SIGN/CD209) and DC-SIGN related (DC-SIGNR/CD299/L-SIGN/CD209L/CLEC4M) are trans-membrane C-type lectins that are calcium-dependent carbohydrate-binding proteins, which can act as cell-adhesion and pathogen-recognition receptors (PMID: 18458800). DC-SIGN and DC-SIGNR share 77% amino acid identity (PMID: 10975799; 11226297). DC-SIGN is expressed on subsets of dendritic cells and macrophages while DC-SIGNR is restricted to endothelial cells of the liver, lymph nodes and placenta (PMID: 17509452). DC-SIGNR serves as a receptor for ICAM-3. It recognizes a wide range of evolutionarily divergent pathogens.

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ untransfected (left) or CD299-transfected HEK-293T (right) cells were surface stained with 0.25 ug Biotin Anti-Human DC-SIGNR/CD299 Rabbit RecAb (Biotin-FcA98205, Clone: 241683G8) or FcZero-rAb™ Biotin Rabbit IgG Isotype Control Recombinant Antibody (Biotin-FcA98136, Clone: 240953C9), and PE Streptavidin (Cat.No. PF00021). Cells were not fixed.

1x10⁶ untransfected (left) or CD299-transfected HEK-293T (right) cells were surface stained with 0.25 ug Biotin Anti-Human DC-SIGNR/CD299 Rabbit RecAb (Biotin-FcA98205, Clone: 241683G8) (red) or FcZero-rAb™ Biotin Rabbit IgG Isotype Control Recombinant Antibody (Biotin-FcA98136, Clone: 240953C9) (blue), and PE Streptavidin (Cat.No. PF00021). Cells were not fixed.