

For Research Use Only

FcZero-rAb™ Biotin Anti-Mouse VE-cadherin/CD144 Rabbit Recombinant Antibody

Catalog Number: Biotin-FcA98251



Basic Information

Catalog Number:

Biotin-FcA98251

Size:

100ug , 500 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0921

GenBank Accession Number:

NM_009868.4

GeneID (NCBI):

12562

UNIPROT ID:

P55284

Full Name:

cadherin 5

Calculated MW:

88kDa

Purification Method:

Protein A purification

CloneNo.:

241953D4

Recommended Dilutions:

FC: 0.25 ug per 10⁶ cells in a 100 µl suspension

Excitation/Emission maxima wavelengths:

-

Applications

Tested Applications:

FC

Species Specificity:

mouse

Positive Controls:

FC : bEnd.3 cells,

Background Information

Cadherins are a family of transmembrane glycoproteins that mediate calcium-dependent cell-cell adhesion and play an important role in the maintenance of normal tissue architecture. Vascular endothelial cadherin (VE-cadherin), also known as Cadherin-5 (CDH5) or CD144, is a member of the type II classical cadherin family of cell adhesion proteins (PMID: 21269602). VE-cadherin is expressed specifically in endothelial cells and mediates homophilic adhesion in the vascular endothelium (PMID: 1522121; 8555485; 21269602). VE-cadherin plays a role in the organization of lateral endothelial junctions and in the control of permeability properties of vascular endothelium (PMID: 1522121). VE-cadherin has also been shown to be required for angiogenesis (PMID: 16473763; 18162609).

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ bEnd.3 cells (left) or C2C12 (right) were surface stained with 0.25 ug Biotin Anti-Mouse VE-cadherin/CD144 Rabbit RecAb (Biotin-FcA98251, Clone: 241953D4) (red) or 0.25 ug Rabbit IgG Isotype Control RecAb (Biotin-FcA98136, Clone: 240953C9) (blue), and PE-conjugated streptavidin. Cells were not fixed.