

For Research Use Only

FcZero-rAb™ Biotin Anti-Human DPP4/CD26 Rabbit Recombinant Antibody

Catalog Number: Biotin-FcA98261



Basic Information

Catalog Number: Biotin-FcA98261	GenBank Accession Number: BC013329	Purification Method: Protein A purification
Size: 100ug , 500 ug/ml	GeneID (NCBI): 1803	CloneNo.: 242045C1
Source: Rabbit	UNIPROT ID: P27487	Recommended Dilutions: FC: 0.25 ug per 10 ⁶ cells in a 100 µl suspension
Isotype: IgG	Full Name: dipeptidyl-peptidase 4	Excitation/Emission maxima wavelengths: -
Immunogen Catalog Number: EG0157	Calculated MW: 88 kDa	

Applications

Tested Applications: FC, ELISA	Positive Controls: FC : human PBMCs,
Species Specificity: human	

Background Information

CD26, also known as DPP4 (dipeptidyl peptidase-4), is a 110 kDa type II cell-surface glycoprotein widely expressed on T cells, activated B cells, activated NK cells and myeloid cells as well as on epithelial cells, fibroblasts, mesothelium, and endothelial cells of a variety of tissues (PMID:33777580; 26919392). It has ecto-enzyme activity in its extracellular domain and cleaves amino-terminal dipeptides with either L-proline or L-alanine at the penultimate position (PMID: 9553764). CD26 plays roles in diverse biological processes such as immunoregulation, glucose homeostasis, and tumorigenesis (PMID: 33777580).



Storage

Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 0.09% sodium azide, pH7.3

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ human PBMCs were surface stained with 0.25 ug Biotin Anti-Human DPP4/CD26 Rabbit RecAb (Biotin-FcA98261, Clone:242045C1) or 5 ul FcZero-rAb™ Biotin Rabbit IgG Isotype Control Recombinant Antibody (Biotin-FcA98136, Clone: 240953C9), and Streptavidin-PE Conjugate (PE-PF00030). Cells were co-stained with FcZero-rAb™ APC Anti-Human CD3 (UCHT1) Rabbit IgG Recombinant Antibody (APC-FcA65570, Clone: UCHT1). Cells were not



1x10⁶ human PBMCs were surface stained with 0.25 ug Biotin Anti-Human DPP4/CD26 Rabbit RecAb (Biotin-FcA98261, Clone:242045C1) and Streptavidin-PE Conjugate (PE-PF00030) (red), or FcZero-rAb™ Biotin Rabbit IgG Isotype Control Recombinant Antibody (Biotin-FcA98136, Clone: 240953C9) (blue). Cells were not fixed.